

Standard Lithiation–Borylation

A user's guide

Table of Contents

Experimental Setup and Handling of Reagents	2
Preparation of Dry Solvents.....	2
Using Acros Seal Bottles	3
General Setup of a Lithiation–Borylation Reaction	5
Lithiation–Borylation of Primary Carbamates/Benzoates.....	7
Lithiation-Borylation of Dialkyl Benzoates	9
Lithiation-Borylation of Secondary Benzylic Carbamates	10
General Procedure: Matteson Homologation	12
General Procedure: Zweifel Olefination	13
Lithiation–Borylation Equivalents Flowchart	14
Reaction Monitoring by ¹¹ B NMR Spectroscopy.....	15
Boronic Ester Handling and Purification.....	16
Volatile Boronic Esters.....	16
Pinacol Boronic Esters	16
Separating Mixtures of Secondary and Tertiary Pinacol Boronic Esters	16
Neopentyl Glycol Boronic Esters	16
Transesterification from tertiary neopentyl glycol to pinacol boronic esters.....	17
Ethylene Glycol Boronic Esters	17
General Procedures for the Preparation of Carbamates and Benzoates	18
Preparation of Primary Carbamates in the Microwave.....	18
Preparation of Primary Carbamates under Reflux – for Big-Scale Reactions.....	18
Selective Carbamate Formation on Primary alcohol using NaH.....	18
Preparation of Primary Benzoates using Phase Transfer Conditions	19
Preparation of Primary Benzoates using the Mitsunobu Reaction	19
Preparation of Secondary Benzylic Carbamates	19
Preparation of Secondary Benzoates	20

Experimental Setup and Handling of Reagents

Preparation of Dry Solvents

In Aggarwal lab, the three most commonly used solvents for lithiation–borylation reactions are: Et₂O, cyclopentyl methyl ether (CPME), and *tert*-butyl methyl ether (TBME). They all are dried using the following procedure. Toluene and CHCl₃ can be dried as well.

1. Take a Young's tube or Strauss flask and add 3Å molecular sieves of an amount corresponding to ¼ of the total volume. Note: you can reuse the molecular sieves as long as you have removed all solvent on high vacuum prior to the introduction of fresh solvent.

2. Heat the open vessel in the microwave for 30 s; wipe away any moisture that has condensed at the neck of the flask; shake the flask to disperse hotspots.

DO NOT put molecular sieves in the microwave if you believe they still contain solvent because this will cause a fire.

The vessel must be open when in the microwave; heating a closed system will cause an **EXPLOSION**.

3. Repeat step 2 until no more water comes out of the sieves as indicated by the amount of moisture condensing at the neck; then place the flask under high vacuum to cool.

4. Once the flask and sieves have cooled (this can take up to an hour) evacuate and backfill the flask three times with N₂. At this point you can use the flask to collect dry solvents from a Grubbs system or fill it with a solvent of your choice directly from the bottle (step 5).

5. Remove the Young's tap from the flask and insert a small funnel into the opening. Bend a needle and insert it next to the funnel to allow gas to escape during addition of solvent.

6. Fill the flask with solvent, remove the funnel and needle, and seal the flask with a Suba Seal septum.

7. Degas your solvent (to remove air introduced while filling the vessel) by vigorously bubbling N₂ through the solution. The traces of O₂ or CO₂ in the solvent interfere with lithiation-borylation by consuming the organolithium reagents, thus reducing the yield. Degassing is not necessary if you get solvents from the Grubbs columns.

8. Close the vessel with a Young's tap under N₂ and leave the solvent over the molecular sieves for 24 hr to dry.

3Å Molecular sieves are excellent for drying a wide range of solvents (D. B. G. Williams, M. Lawton, *JOC*, **2010**, 75, 8351).

Using Acros Seal Bottles

Parafilm

When wrapping the bottle, ensure that the parafilm covers both the red lid and the joint of the seal to the bottle (see right) – this prevents moisture ingress. Don't just wrap the red cap (left).



Puncturing

As you can see on the left, the seal has been punctured a few times towards the centre; the seal is not yet compromised. However, the example on the right shows how repeated puncturing in similar places towards the centre of the seal seriously compromises its integrity. The sBuLi in the right bottle has gone off, even though is still ~80% full.




How can this be avoided?

- Only puncture the Acros seal around the edge of the seal against the plastic rim. This is the strongest part of the seal.

- Do not puncture the same hole more than once, or too close to an existing one.
- Only use the long green needles (Gauge 21) to dispense reagents from these bottles; the yellow needles are too wide.
- Only puncture once for each portion of reagent. Use your reaction flask or a separate nitrogen-flushed vessel to flush the syringe with inert gas.
- Ensure the reagent bottle is stored upright.
- Make sure the reagent has warmed to room temperature before removing parafilm (20-30 min depending on the bottle size).
- Use a nitrogen-filled balloon and **not your manifold** to keep the bottle under nitrogen. The use of a manifold for such an operation is **VERY DANGEROUS** – if the syringe and needle come apart the reagent will be pushed out until its level is below the needle tip.

Below is an example of good (top-right quadrant) and bad (bottom-left quadrant) ways to puncture an Acros seal.



General Setup of a Lithiation–Borylation Reaction

1. Take a clean Schlenk flask/tube with a stirrer bar and heat it on a Bunsen burner until the flame becomes slightly orange. Avoiding heating the side-arm tap too much as it may crack. Insert a Suba-Seal septum and allow the flask to cool under vacuum.

- The stirrer must stir freely at the bottom of the Schlenk flask and not be wider than the diameter of the vessel

- Using a gas adapter to dry your Schlenk is unnecessary – Suba-Seal septa are cheap and work just fine.

- Parafilming the Suba-Seal septum is not necessary – if you think it leaks, better throw it away and take a new one.

2. Once the flask has cooled evacuate and backfill the flask with N_2 three times.

- Even though N_2 supplies are typically dry I would advise having two drying traps attached to your N_2 supply – one with molecular sieves for coarse drying and the other with Drierite with some self-indicating silica for fine drying.

3. For carbamates or benzoates, which are oils, fill the dead volume of a syringe with your reagent and tare the syringe on the balance. After the balance has settled draw up more reagent and remove any bubbles in the syringe until you have the desired weight. Inject the syringe into the reaction flask, but do not pump the syringe, then return the dead volume to your flask of starting material.

4. For solid benzoates or carbamates, weigh out the desired amount and put the solid into the Schlenk flask; then evacuate the Schlenk and refill with N_2 .

- for steps 3 and 4, minimize the amount of time your Schlenk flask is exposed to the outside atmosphere

5. At this point add your diamine by syringe, TMEDA, (+)- or (–)-sparteine, by syringe using the known densities of these compounds to measure out the correct volumes under N_2 . $d = 0.775$ (TMEDA), 1.02 (sparteine)

Handling of (+)- or (–)-sparteine

- If the sparteine has discoloured to yellow (or sometimes red), redistill it before use – it should be colourless.

- Sparteine should be stored in a Schlenk flask with a glass stopper in the freezer. Flasks with Teflon taps should not be used as the joint can lose its seal at low temperatures.

- Sparteine solidifies in freezer. It should be room temperature before use. You can accelerate the melting by partially immersing the sealed sparteine contained into a beaker with warm water.

- Evacuate and backfill the head space of the tap with dry N_2 /Ar three times, then open the tap.

- Remove the glass stopper and replace it with a septum. Place the whole flask under vacuum, then refill with N_2 .

- Evacuate a syringe three times in a separate vessel with N_2 , then move the syringe to the sparteine vessel and draw up the desired amount. Pull a small N_2 bubble between the sparteine and the needle and transfer the syringe to the reaction vessel, pushing out the bubble first. Pull up a head of N_2 for the remaining dead volume and pump the syringe back into the sparteine vessel.

- Remove the Suba-Seal septum and replace it with a clean glass stopper. Put the sparteine flask under vacuum and refill with N₂. Close the glass tap and parafilm both the end of the tap and the glass stopper joint, then place the flask back in the freezer.

6. Evacuate and backfill a syringe in a separate nitrogen-flushed flask. Evacuate and backfill your Strauss flask containing your reaction solvent and replace the Young's tap with a Suba-Seal septum. Draw up your solvent, pull up a cushion of N₂ and transfer the solvent to the reaction flask, washing the sides of the vessel to dislodge reagents that have adhered to the sides.

7. Ensure your carbamate/benzoate and diamine are completely dissolved. Place the reaction flask in a cooling bath and ensure that the stirrer moves freely without splashing the reaction mixture all over the inside of the Schlenk flask.

8. Fill the bath with dry ice and then acetone and allow the flask to cool for ~2 mins. The bath should be full of dry ice; a bit of CO₂ at the bottom will not suffice.

9. Attach a nitrogen-filled balloon to a bottle of sBuLi. Fill a syringe with nitrogen in a separate N₂-filled flask, then insert it into the BuLi bottle. Draw up the required amount of sBuLi. Rotate the syringe upside-down (with the needle still inside the BuLi bottle), pull up a protective cushion of N₂, then carefully move the syringe to your reaction vessel.

- Do not evacuate the syringe in the sBuLi bottle as this shortens the seal lifetime.

- Always pierce the Acros seal around the edge of the seal and not in the middle.

10. Add the sBuLi dropwise (~1 drop every 1-2 seconds) until the addition is complete

- Push out the head of N₂ (with the syringe upside down) before adding sBuLi itself, so that you do not accidentally add the dead volume.

- Make sure that sBuLi drops directly into the reaction solution: if it hits the wall of the vessel, it will freeze

-----Wait the required lithiation time-----

(2 - 5 h for primary carbamates/benzoates, 0.25 - 1 h for secondary benzylic/allylic/propargylic carbamates, 2-8 h for dialkyl benzoates)

11. If the boronic ester is an oil, weigh it out as per step 3, and add the boronic ester neat in a dropwise manner.

- Do not add the dead volume of the needle, return this to your flask of reagent.


- If the ice bath is bubbling vigorously, you are adding the boronic ester too quickly – lots of heat is generated.

12. If the boronic ester is a solid, weigh it out as per step 4 into a flame-dried vial under N₂, and dissolve in the reaction solvent to 1M concentration. Then add this solution dropwise to the reaction mixture.

- If the ice bath is bubbling vigorously, you are adding the boronic ester too quickly – lots of heat is generated.

13. Give enough time for the formation of the boronate complex (typically < 1h, however more hindered substrates may require up to 3h), then follow the procedure appropriate for your particular carbamate/boronic ester combination to effect 1,2-migration.

Lithiation–Borylation of Primary Carbamates/Benzoates



Appropriate citations for this methodology

ACIE, **2007**, 46, 7491; *Chem. Commun.*, **2011**, 47, 12592; *Org. Synth.*, **2011**, 88, 247; *Acc. Chem. Res.*, **2014**, 47, 3174.

Carbenoid as Limiting Reagent

To a solution of carbamate/benzoate (1.0 mmol) and (–)-sparteine (1.3 mmol) in dry Et₂O (5 mL) at –78 °C was added sBuLi (1.3 M in cyclohexane, 1.3 mmol) dropwise. This mixture was stirred at –78 °C for the appropriate lithiation time (most substrates will completely lithiate within 5 hours). A solution of boronic ester (1 M in Et₂O, 1.3 mmol) was added dropwise and the reaction mixture was stirred for a further hour at –78 °C. The reaction mixture was then heated at reflux until the boronate complex was consumed. The reaction mixture was diluted with Et₂O (20 mL) and quenched through addition of 2M aq. HCl (20 mL). After stirring for 10 min, the phases were separated, and the aqueous phase was re-extracted with diethyl ether (2 × 20 mL). The combined organic phases were dried over anhydrous MgSO₄, filtered, and the solvent was removed in vacuo. The aqueous phase was retained for recovery of sparteine [see - *Org. Synth.* **1997**, 74, 23, *Nature*, **2014**, 513, 183]. The crude product was purified by column chromatography to afford the secondary boronic ester.

To isolate the secondary alcohol:

The reaction mixture was diluted with THF (5 mL) and then cooled to 0 °C. A premixed solution of NaOH (2 M)/H₂O₂ (30%) (2:1 v/v, 4.0 mL) was added dropwise. The biphasic mixture was stirred for 1h at rt and quenched with 2 M NaOH (10 mL). The layers were separated and the aqueous layer was extracted with Et₂O (3 x 15 mL). The combined organic layers were dried (MgSO₄) and concentrated *in vacuo*. The crude product was purified by flash chromatography.

Formation of MgBr₂•Et₂O

In some cases, addition of this Lewis acid is necessary to promote 1,2-migration with primary carbamates.

1,2-Dibromoethane (1.3 mmol, 1.3 eq.) was added to Mg turnings (2 eq., 2.0 mmol) suspended in Et₂O (5 mL). The reaction mixture was gently heated to initiate the reaction and stirred until reflux has stopped. The biphasic mixture was added to the lithiation-borylation reaction mixture, stirred at –78 °C for 20 min, then warmed to r.t. and heated at reflux for ≥ 12h.

Boronic Ester as Limiting Reagent


To a solution of carbamate/benzoate (1.5 mmol) and (–)-sparteine (304 mg, 1.5 mmol) in dry Et₂O (7.5 mL) at –78 °C was added sBuLi (1.3 M in cyclohexane, 1.4 mmol) dropwise. This mixture was stirred at –78 °C for the appropriate lithiation time. A solution of boronic ester (1 M in Et₂O, 1.0 mmol) was added dropwise and the reaction mixture was stirred for a further hour at –78 °C. The reaction mixture was then heated at reflux until the boronate complex was consumed. The reaction mixture was diluted with Et₂O (20 mL) and quenched through addition of 2M aq. HCl (20 mL). After stirring for 10 min, the phases were separated, and the aqueous phase was re-extracted with diethyl ether (2 × 20 mL). The combined organic phases were

dried over anhydrous MgSO_4 , filtered, and the solvent was removed in vacuo. The aqueous phase was retained for recovery of sparteine. The crude product was purified by column chromatography to afford the secondary boronic ester.

Notes:

- The lithiated carbamate should appear pale yellow and the benzoate should appear pale yellow, brown or purple depending on how it has been prepared.
- As benzoates and sparteine are quite viscous it is important to ensure they are fully dissolved before introduction of the organolithium.

Lithiation-Borylation of Dialkyl Benzoates



Appropriate citations for this chemistry

JACS, **2013**, 135, 16054; Acc. Chem. Res., **2014**, 47, 3174.

To a vigorously stirred solution (without splashing) of benzoate (0.50 mmol) and TMEDA

(0.46 ml, 3.00 mmol) in anhydrous CPME (3 ml) at $-60\text{ }^{\circ}\text{C}$ (internal temperature) under a nitrogen atmosphere, was added sBuLi (1.3 M in hexane, 0.62 ml, 0.80 mmol) dropwise over 10 min. After 2 h, a solution of the neopentylglycol boronic ester (1.0 mmol) in CPME (0.5 ml) was added dropwise over 10 min. The reaction mixture was stirred at $-60\text{ }^{\circ}\text{C}$ for 1 h and then the cooling bath was removed and the reaction mixture was stirred at the $50\text{ }^{\circ}\text{C}$ (migration temperature) overnight (~ 16 h).

For isolation of the tertiary pinacol boronic ester -

The volatiles were removed from the reaction mixture in vacuo and the residue dissolved in THF (6 ml). To the stirred solution was added pinacol (709 mg, 6.00 mmol) and 0.2 M aq. NaOH (2 ml). The reaction mixture was stirred for 5 h at $60\text{ }^{\circ}\text{C}$ at which point GCMS showed 95% conversion of the tertiary neopentylglycol boronic ester to the corresponding pinacol boronic ester. The reaction mixture was cooled to ambient temperature and Et_2O (10 ml) and water (5 ml) was added. The phases were separated and the aqueous phase extracted with Et_2O (3×20 ml). The combined organic phases were washed with brine, dried over MgSO_4 and the solvent removed in vacuo. The crude material was purified by flash column chromatography eluting 1–1.8% Et_2O /pentane to give tertiary boronic ester.


For isolation of the tertiary alcohol -

The reaction mixture was cooled to $0\text{ }^{\circ}\text{C}$ (water/ice bath) and degassed THF (2 ml) containing BHT (~ 2 mg) was added. An ice-cold degassed mixture of 3 M NaOH (3.6 ml) containing EDTA (1 g/L) and 30% aqueous H_2O_2 (1.8 ml) was added all at once. The reaction mixture was vigorously stirred and allowed to reach ambient temperature overnight. The reaction mixture was diluted with water (5 ml) and extracted with Et_2O (4×20 ml). The combined organic phases were washed with brine and dried over MgSO_4 . The solvent was removed in vacuo and the crude material purified by flash chromatography.

Notes:

- It is important the reaction mixture does not splash all over the flask during the lithiation step, otherwise the lithiation is not as successful
- $-60\text{ }^{\circ}\text{C}$ should be maintained using a pre-equilibrated cryostat bath using acetone as the coolant. A mechanical stirrer should be used to ensure the cooling bath is homogeneous, DO NOT use a stirrer bar to stir the bath as this will cause the reaction flask to stir erratically.
- Attempting to put two reaction vessels into the $-60\text{ }^{\circ}\text{C}$ cooling bath at once will lead to less reproducible results due to poor stirring.
- Exact lithiation times are substrate specific see - J. Am. Chem. Soc., **2013**, 135, 16054
- The lithiated benzoate should appear purple or brown.
- The tertiary pinacol boronic ester products and the starting material secondary benzoates will be very similar in polarity, unless the boronic ester used contains polar groups.

Lithiation-Borylation of Secondary Benzylic Carbamates



Appropriate citation for this chemistry

Nature, **2008**, 456, 778; *ACIE*, **2012**, 49, 5142; *Chem. Sci.*, **2015**, 6, 3718.

A solution of secondary benzylic carbamate (1.00 mmol, 1.0 equiv) in anhydrous diethyl ether or TBME (3.0 mL) was cooled to -78 °C. sBuLi (1.00 mL, 1.30 mmol, 1.3 equiv) was added dropwise and the reaction mixture was stirred at this temperature for 0.15–1 h. A solution of the boronic ester (1.50 mmol, 1.5 equiv) in anhydrous diethyl ether (1.5 mL) was added dropwise and the mixture was stirred for 1–3 h at -78 °C [boronate complex formation]. A 1.0 M solution of MgBr₂ in anhydrous MeOH (1.30 mL, 1.30 mmol, 1.3 equiv) was added slowly at -78 °C. After 5 min, the cooling bath was removed and stirring was continued at room temperature until disappearance of the boronate complex (5–8 ppm) monitored by ¹¹B NMR spectroscopy. The reaction mixture was then cooled to 0 °C and 1.0 M aq. KH₂PO₄ (2.0 mL) was added slowly. After stirring for 10 min, the phases were separated, and the aqueous phase was re-extracted with diethyl ether (4 × 20 mL). The combined organic phases were dried over anhydrous MgSO₄, filtered, and the solvent was removed *in vacuo*. The crude product was purified by column chromatography to afford the tertiary boronic ester.

For isolation of the tertiary alcohol -

The reaction was cooled to 0 °C (water/ice bath) and degassed THF (2 mL) containing BHT (~2 mg) was added. An ice-cold degassed mixture of 3 M NaOH (3.6 mL) containing EDTA (1 g/L) and 30% aqueous H₂O₂ (1.8 mL) was added all at once. The reaction mixture was vigorously stirred and allowed to reach ambient temperature overnight. The reaction was diluted with water (5 mL) and extracted with Et₂O (4 × 20 mL). The combined organic phases were washed with brine and dried over MgSO₄. The solvent was removed *in vacuo* and the crude material purified by flash chromatography.

Notes:

- For carbamates with electron-donating groups, TMEDA (1.3 eq) should be added before performing your lithiation.
- The lithiated carbamate should appear anywhere from yellow to deep orange in colour.
- THF should not be used for these lithiations as it will cause your lithiated carbamate to become configurationally unstable.
- MgBr₂ in MeOH was prepared by drying MgBr₂ at 100 °C overnight under high vacuum, followed by addition of anhydrous MeOH. A large batch can be prepared and stored over a long period of time; however, take care as during the initial stages of dissolution, the MeOH will start to boil. Any solids formed (MgO) should be left to settle before using the reagent.


For Pyridyl Boronic Ester Substrates

Carbamate (0.40 mmol) was dissolved in anhydrous Et₂O (1 ml), cooled to -78 °C and sBuLi (0.37 ml of a 1.3 M hexane:cyclohexane solution, 0.48 mmol) was added dropwise over 5 min. After stirring for 20 min, the respective pyridyl boronic ester (0.60 ml of a 1.0 M Et₂O solution, or 1.20 ml of a 0.5 M PhMe solution 0.60 mmol) was added dropwise over 5 min. The reaction mixture was stirred for an additional 1 h, then MgBr₂ (0.60 ml of a 1.0 M MeOH solution) was added dropwise over 2-3 min. After 15 min the cooling bath was removed and stirring was continued at rt for 16 h. The reaction mixture was diluted with anhydrous THF (2 ml), chilled to ~0 °C, and 0.30 ml of a degassed mixture of 3.0 M aq. NaOH (containing EDTA, 1.0 g/l) and 30% aq. H₂O₂, 1:1 (vol) was added. The cooling bath was removed and the reaction mixture was stirred at rt for 30 min. The reaction was diluted with water (5 ml) and Et₂O (5 ml) and extracted with Et₂O (3 × 5 ml). The combined organic phases were washed with water (5 ml), brine (10 ml), dried over MgSO₄, filtered and concentrated *in vacuo*. The crude product was purified by flash chromatography to give pure tertiary alcohol.

For Homologation of Tertiary Boronic Esters - *J. Am. Chem. Soc.*, **2014**, *136*, 17370

Tertiary boronic ester (0.3 mmol) was dissolved in anhydrous Et₂O (1 ml) and cooled to 0 °C. MeMgBr (0.6 ml of a 1 M Et₂O solution, 0.60 mmol) was added and the reaction allowed to stir at rt for 5 min. Secondary benzylic carbamate (0.2 mmol) was dissolved in anhydrous Et₂O (1 ml), cooled to -78 °C and sBuLi (160 µl of a 1.30 M hexane:cyclohexane solution, 0.21 mmol) added dropwise over 2 min. After stirring for 20 min, both the lithiated carbamate and borane were cooled to -96 °C (MeOH/N₂(l) bath) and the borane transferred into to the lithiated carbamate via cannula. The reaction mixture was stirred at this temperature for 5 min, then rt for 1 h. The reaction mixture was chilled at 0 °C, and 0.30 ml of a degassed mixture of 3.0 M aq. NaOH (containing EDTA, 1.0 g/l) and 30% aq. H₂O₂, 1:1 (vol) was added. The cooling bath was removed and the reaction mixture was stirred at rt for 30 min. The reaction was diluted with water (5 ml) and Et₂O (5 ml) and extracted with Et₂O (3 × 5 ml). The combined organic phases were washed with water (5 ml), brine (10 ml), dried over MgSO₄, filtered and concentrated *in vacuo*. The crude product was purified by flash column chromatography (EtOAc/petroleum ether 40-60, 1:9) to give pure tertiary alcohol.


General Procedure: Matteson Homologation



Appropriate citations

JACS, **2015**, 137, 4398; *ACIE*, **2011**, 50, 3760; *JACS*, **1980**, 102, 7588; *Organometallics*, **1983**, 2, 1529, *Organometallics*, **1992**, 11, 1948.

In a flamed-dried flask under nitrogen a solution of boronic ester (1.0 equiv.) and chloriodomethane (3.0 equiv.) in anhydrous Et₂O (0.2 M) was cooled to –95 °C (MeOH/liq.N₂). *n*-Butyllithium (1.5–1.6 M in hexanes, 2.95 equiv.) was added slowly using a syringe pump (around 3 drops per minute). The reaction mixture was stirred for 30 min at –95 °C and then warmed to RT and stirred for a further 1 h. The cloudy solution was filtered through a plug of silica (~ 10 mm depth of silica, using a glass-fritted Büchner funnel connected directly to a receiving round-bottom flask) to give a colourless solution. The silica was washed with Et₂O (reagent grade), the funnel was removed and solvent was evaporated under reduced pressure to give the crude boronic ester (for the example shown above: 99%, <1% starting material by GCMS analysis).




In a flamed-dried flask under nitrogen a solution of boronic ester (1.0 equiv.) and bromochloromethane (3.0 equiv.) in anhydrous Et₂O (0.25 M) was cooled to –78 °C (acetone/dry ice). *n*-Butyllithium (1.5–1.6 M in hexanes, 2.5 equiv.) was added slowly using a syringe pump (around 3 drops per minute). The reaction mixture was stirred for 30 min at –78 °C then warmed to RT and stirred for a further 1 h. The cloudy solution was filtered through a plug of silica (~ 10 mm depth of silica, using a glass-fritted Büchner funnel connected directly to a receiving round-bottom flask) to give a colourless solution. The silica was washed with Et₂O (reagent grade), the funnel was removed and solvent was evaporated under reduced pressure to give the crude boronic ester (for the example shown above: 99%, <1% starting material by GCMS analysis).

General Procedure: Zweifel Olefination

Appropriate citations


ACIE, **2011**, 50, 3760; *ACIE*, **2009**, 48, 6317; *JACS*, **1967**, 89, 3652; *JOC*, **1976**, 41, 3947.

Use a vinyl Grignard for simple unhindered substrates:



To a solution of boronic ester (0.30 mmol) in anhydrous THF (3.0 mL, 0.1 M) at room temperature was added a vinylmagnesium bromide (1.0 M in THF, 1.2 mmol) dropwise. The resulting mixture was stirred at room temperature for 30 min and cooled down to -78°C . A solution of iodine (1.2 mmol) in MeOH (4 mL) was added dropwise to the reaction mixture via cannula; 30 min later, a solution of MeONa (2.4 mmol) in MeOH (5 mL) was added. The reaction mixture was then allowed to warm to room temperature and stirred for an additional 1 h, diluted with pentane (40 mL) and washed with a 20% aq. $\text{Na}_2\text{S}_2\text{O}_3$ solution (5 mL) and water (10 mL). The phases were separated, the aqueous layer was extracted with pentane (2×20 mL); the combined organic layers were washed with brine (10 mL), dried (MgSO_4), filtered and concentrated *in vacuo*. The crude product was purified by flash column chromatography.

Use vinyl lithium for more hindered substrates:



To neat tetravinyltin (0.30 mmol) at room temperature was added *n*BuLi (1.6 M in hexanes, 0.60 mmol) dropwise. The resulting solution was stirred for 5 min at room temperature, and then allowed to settle for 3 h. The supernatant was removed, and the white solid was carefully washed with hexane (2 \times), attention being paid not to remove any of the precipitate. The solid vinyl lithium was then dissolved in THF (0.6 mL, 1 M, 0.6 mmol).

To a solution of boronic ester (0.30 mmol) in anhydrous THF (3.0 mL, 0.1 M) at -78°C was added dropwise the vinyl lithium solution. The resulting mixture was stirred for 45 min, removed from the bath and stirred at room temperature for 30 min. The reaction mixture was then cooled down to -78°C and a solution of iodine (0.6 mmol) in MeOH (2 mL) was added dropwise to the reaction mixture; 30 min later, a solution of MeONa (2.4 mmol) in MeOH (5 mL) was added. The reaction mixture was stirred at -78°C for 30 min, then allowed to warm to room temperature and stirred for an additional 1 h, then diluted with pentane (40 mL) and washed with a 20% aq. $\text{Na}_2\text{S}_2\text{O}_3$ solution (5 mL) and water (10 mL). The phases were separated and the aqueous layer was extracted with pentane (2×20 mL); the combined organic layers were washed with brine (10 mL), dried (MgSO_4), filtered and concentrated *in vacuo*. The crude material was purified by flash column chromatography to yield the desired product.

Notes:

The most important stage of the reaction is the warming of the reaction mixture to room temperature after addition of the vinyl lithium. This ensures that you achieve complete boronate complex formation.


We have reported alternate procedures involving the addition of just I_2 in MeOH to effect 1,2-metallated rearrangement and elimination. However, investigation of more complex and challenging substrates has indicated that the addition of NaOMe is beneficial and therefore should be used more generally.

Lithiation–Borylation Equivalents Flowchart

Deciding on the appropriate relative number of equivalents of substrates and reagents for your lithiation–borylation reaction can sometimes be difficult:

You should first identify your more precious reagent, the boronic ester or the lithiated carbamate/benzoate; this consideration is most important for total synthesis as it will determine your longest linear sequence.



If you are preparing starting materials for methodology, your boronic ester will be your most precious compound. You will want to ensure it is all consumed so you don't have to separate homologated boronic ester from starting boronic ester, a process that can be difficult. The next decision is whether to use a benzoate and carbamate as your carbenoid precursor. If you expect that your product boronic ester will be very apolar then you should use a carbamate as it will be more polar than your product. However, if your boronic ester is going to be polar then use the benzoate as it will typically be more apolar than your product. This analysis assumes that the 1,2-migration step in your lithiation-borylation reaction is facile. When 1,2-migration is challenging (phenyl, vinyl, Me, EWGs \rightarrow to the boronic ester are all poor migrating groups), the use of the benzoate gives superior results.



Reaction Monitoring by ^{11}B NMR Spectroscopy

Lithiation–borylation reactions can be monitored by ^{11}B NMR spectroscopy. Your starting material boronic ester has a chemical shift of ~ 30 ppm; when you form a boronate complex, the chemical shift changes to around 0 to 10 ppm. After 1,2-migration, the product boronic ester has a chemical shift of ~ 30 ppm, thus allowing you to monitor the consumption of your boronate complex by ^{11}B NMR spectroscopy. On its own, ^{11}B NMR data do not tell you if your intermediate boronate complex has reverted to starting materials or has undergone 1,2-migration to products (in both cases, a boronic ester is formed); however, it does inform you as to whether your boronate complex has been consumed, and in the case where it has, your reaction is complete (for better or for worse).

GCMS and TLC analysis of your crude reaction mixture can be used in tandem with your ^{11}B NMR analysis to determine if starting material or product has formed.




A useful resource for other chemical shifts can be found here

<http://www.chemistry.sdsu.edu/research/BNMR/>

Boronic Ester Handling and Purification

Preparation of Simple Boronic Ester Starting Materials from the Corresponding Boronic Acids



Boronic acid (1.0 mmol) and diol (1.0 mmol) were stirred in anhydrous Et₂O (1.0 ml) at ambient temperature under a nitrogen atmosphere for 2–10 h. Flame-dried MgSO₄ (5 mmol) was added to the reaction mixture and stirring continued for 12 h. The reaction mixture was filtered and the solvent removed in vacuo. The crude material was purified by distillation or flash column chromatography.

Volatile Boronic Esters

If your product boronic ester is of low molecular weight, take care when evaporating your solvent as the boronic ester can azeotrope with the Et₂O. The best way to avoid this is to evaporate the Et₂O without using the vacuum on your rotary evaporator and setting the bath to 40–50 °C – this works best with an efficient dry-ice condenser.

Pinacol Boronic Esters

Pinacol boronic esters are, for the most part, bench stable and will only slowly decompose through oxidation/protodeboronation over periods of time >1 year. They can be purified by column chromatography or Kugelrohr distillation. The latter is especially useful when your starting material boronic ester is substantially more volatile than your product boronic ester. Do not be afraid to heat your boronic esters in the Kugelrohr apparatus; many are stable to such treatment (See the Supporting Information of the following paper: *Angew. Chem. Int. Ed.*, **2010**, 49, 5142).

Some substrates require extra care in their handling such as allylic or propargylic boronic esters as they can decompose upon chromatographic purification; in some cases, switching the constituent polar eluent from diethyl ether to ethyl acetate (or indeed vice versa) can lead to superior yields of product.

Separating Mixtures of Secondary and Tertiary Pinacol Boronic Esters

In a few isolated cases where tertiary pinacol boronic esters are contaminated with small amounts of starting secondary pinacol boronic ester, the latter could be removed by the following treatment: the mixture was dissolved in diethyl ether; the organic phase was washed several times with 0.5 M aq. NaOH solution; the solvent was removed in vacuo to give the pure tertiary pinacol boronic ester.

Note: This was done with specific hindered tertiary boronic esters and it is not known how general this procedure is.

Neopentyl Glycol Boronic Esters

These boronic esters must be stored under N₂ if they are to be kept for long periods of time, otherwise oxidation/protodeboronation will occur. For instance, a good way to store simple boronic esters, such as *i*PrBneo, is in a Young's tube.

Neopentyl glycol boronic esters can be purified by chromatography, but the amount of time the material is in contact with the silica gel needs to be minimal. By TLC, the spots will usually appear as streaks for unhindered neopentyl glycol boronic esters and well-defined spots for hindered neopentyl glycol boronic esters. Usually, simple primary neopentyl glycol esters will decompose on silica gel faster than more

hindered secondary/tertiary neopentyl glycol boronic esters, thus allowing hindered boronic ester products to be separated from less hindered starting boronic esters through the latter's more rapid decomposition on silica gel.

Transesterification from tertiary neopentyl glycol to pinacol boronic esters

After a lithiation–borylation reaction (2.00 mmol scale), the volatiles were removed from the reaction mixture *in vacuo* and the residue dissolved in THF (6 ml). To the stirred solution was added pinacol (709 mg, 6.00 mmol) and 0.2 M aq. NaOH (2 ml). The reaction mixture was stirred for 5 h at 60 °C or until such point as the tertiary neopentyl boronic ester is converted into the corresponding pinacol boronic ester (as determined by GCMS or otherwise). The reaction mixture was cooled to ambient temperature and Et₂O (10 ml) and water (5 ml) was added. The phases were separated and the aqueous phase extracted with Et₂O (3 × 20 ml). The combined organic phases were washed with brine, dried over MgSO₄ and the solvent removed *in vacuo*. The crude material was purified by flash column chromatography.

Ethylene Glycol Boronic Esters


These boronic esters are exceptionally air and moisture sensitive and should always be stored under N₂. They should be distilled or used without purification as they cannot be purified by column chromatography. These boronic esters are very hygroscopic so solid ethylene glycol boronic esters should be stored as solutions over CaCl₂ or 3 Å MS, and oils should be stored neat over CaCl₂.

Transesterification from ethylene glycol to pinacol boronic esters after a lithiation–borylation reaction

After a lithiation–borylation reaction, the volatiles were removed *in vacuo* (~0.3 mbar) and the residue was redissolved in anhydrous Et₂O. Pinacol (1.0–1.5 eq.) was added and the mixture was stirred at room temperature for 16 hours before being diluted with H₂O. The phases were separated and the aqueous phase extracted with Et₂O. The combined organic phases were washed with brine, dried over MgSO₄ and the solvent removed *in vacuo*. The crude material was purified by flash column chromatography.

General Procedures for the Preparation of Carbamates and Benzoates

Preparation of Primary Carbamates in the Microwave




N,N-diisopropylcarbamoyl chloride (1.2 eq.), triethylamine (1.3 eq.) and primary alcohol (1 eq.) were added to a flame-dried microwave vial and dissolved in toluene (1M in alcohol). The mixture was stirred briefly to ensure that everything was

dissolved and then heated in the microwave reactor for 1–2 h at 150 °C. After cooling to room temperature, the reaction mixture was filtered through a short plug of silica to remove Et₃N•HCl and the filter cake washed with Et₂O. Further purification of the carbamate by standard column chromatography will afford analytically pure material.

Notes:


- Most substrates will be transformed within 1 h at 150 °C. Sometimes conversion is incomplete after such time owing to the precipitated Et₃N•HCl preventing efficient mixing; in such cases, shaking the vial and heating the mixture again for 1 h at 150 °C should suffice to achieve complete conversion.
- For substrates which fail to lithiate to completion, additional distillation of the carbamate after chromatography may be required and/or azeotropic removal of trace water using dry toluene.
- The biggest vials that fit in our microwave can take 20 mL which equates to roughly 15 mL of toluene + other reagents, and is therefore ideally suited to small-to-medium-scale reactions.

Preparation of Primary Carbamates under Reflux – for Big-Scale Reactions



N,N-diisopropylcarbamoyl chloride (1.2 eq.), triethylamine (1.3 eq.) and primary alcohol (1 eq.) were added to a flame-dried round-bottomed flask and dissolved in dichloromethane (1M in alcohol). A reflux condenser was attached and the reaction mixture was heated at reflux for 24 to 48 h. After cooling to

room temperature, the reaction mixture was quenched through the addition of water and the phases were separated. After extracting the aqueous phase with DCM (x 3), the combined organics were dried over MgSO₄, filtered and concentrated. Purification of the residue by column chromatography or in some cases kugelrohr distillation will afford analytically pure material.




Selective Carbamate Formation on Primary alcohol using NaH

(*S*)-(+)-1,3-Butane diol (3.0 mL, 33.5 mmol, 1 eq.) was added dropwise to a suspension of NaH (60% dispersion in mineral oil, 2.74 g, 68.6 mmol) in THF (330 mL) at 0 °C. The mixture was warmed to RT and stirred for 2 h. The mixture was cooled to 0 °C and a solution of *N,N*-diisopropylcarbamoyl chloride (5.48 g, 33.5 mmol) in THF (33 mL) was added dropwise. The mixture was allowed to warm to RT slowly over 16 h. Water (300 mL) and Et₂O (300 mL) were added and the layers were separated; the aqueous phase was extracted with Et₂O (2 × 300 mL). The combined organic phases were washed with brine (300 mL), dried (MgSO₄) and concentrated in vacuo. Purification by column chromatography (pentane/EtOAc 3:1) gave (*S*)-3-hydroxybutyl diisopropylcarbamate as an oil (5.41 g, 74%).

Notes: - Can also be performed on simple primary alcohols however the CbCl/Et₃N is practically simpler for these substrates

- Add the alcohol slowly otherwise it will foam out of the flask

Preparation of Primary Benzoates using Phase Transfer Conditions




A biphasic mixture of 2,4,6-triisopropylbenzoic acid (20.2 g, 81.3 mmol, 1.0 eq), $\text{NBu}_4(\text{HSO}_4)$ (2.21 g, 6.5 mmol, 0.08 eq), NaOH (10.1 g, 252.0 mmol, 3.1 eq) and bromoethane (30.0 mL, 407 mmol, 5.0 eq) in CHCl_3 (400 mL) and H_2O (320 mL) in a 1 L round-bottomed flask was vigorously

stirred overnight at room temperature. The phases were separated and the aqueous phase extracted with CH_2Cl_2 (3 × 100 mL). The combined organic phases were washed with brine (300 mL), dried (MgSO_4) and concentrated under reduced pressure. The crude product was dissolved in pentane (60 mL) and the insoluble salts filtered off. The solvent was removed from the filtrate under reduced pressure to give ethyl TIB ester (20.1 g, 90%) as a colourless oil.

Notes:

- Works for simple aliphatic alkyl halides such as ethylbromide, *n*-propyl bromide, 3-butenyl bromide, 1-bromo-3-phenyl-propane, but not for more hindered substrates such as isobutyl bromide.
- Can be used without further purification.

Preparation of Primary Benzoates using the Mitsunobu Reaction




To a stirred solution of PPh_3 (11 mmol), secondary alcohol (10 mmol) and 2,4,6-triisopropylbenzoic acid (11.5 mmol) in THF (15 mL) at 0 °C (ice bath), was added DIAD (11 mmol) dropwise over 10 min. After stirring the mixture for 4 h at 0 °C (ice bath), the volatiles were removed *in vacuo*.

The residue was dissolved in pentane (15 mL) and the resulting solution stirred in for 5 min. The white suspension was filtered and the filter cake washed with pentane (100 mL). The solvent was removed *in vacuo* and the residue purified by flash column chromatography to give pure benzoate.

Notes:

- Separation of trace PPh_3 can be difficult sometimes; however, the purification can be simplified through dissolving the triturated crude mixture in THF, adding an excess of 30% aq. H_2O_2 , and stirring the resulting mixture, thus converting PPh_3 into OPPh_3 .
- Developed TLC plates will often suggest an easy purification by column chromatography; however, the use of too polar an eluent will cause OPPh_3 to slowly leach through the column and contaminate product fractions.

Preparation of Secondary Benzylic Carbamates



Reduction of benzylic ketones

To a solution of ketone (10.0 mmol, 1 eq.) in MeOH (7 mL) and THF (10 mL) at 0 °C was added NaBH_4 (567 mg, 15.0 mmol, 1.5 eq.) portionwise over 5 minutes with vigorous stirring. The reaction mixture was then warmed to RT and stirred for 30 min at which point TLC (20% EtOAc:petrol) indicated complete conversion

of starting material. The reaction was quenched by addition of aq. NH_4Cl (5 mL) and diluted with H_2O (10 mL) and EtOAc (100 mL). The layers were separated and the organic layer was washed sequentially with H_2O (15 mL) and brine (2×15 mL). The organic layer was dried over MgSO_4 , filtered and concentrated *in vacuo* to give essentially pure product.


Enzymatic resolution of benzylic alcohols

To a solution of benzylic alcohol (9.8 mmol, 1 eq.) in diisopropyl ether (4 mL) was added acrylic resin bound lipase from *Candida Antarctica* (59 mg, 6 mg per mmol of alcohol) followed by vinyl acetate (4.3 mL, 49 mmol, 5 eq.). The suspension was then heated to 50 °C and stirred for 16 h, after which chiral HPLC analysis should indicate the alcohol having >99:1 *er*. The reaction mixture was filtered through a plug of SiO_2 with EtOAc as eluent, concentrated *in vacuo* and purified by flash column chromatography (20% EtOAc:petrol) to give the enantioenriched (*S*)-alcohol and (*R*)-acetate products.

Carbamoylation of secondary benzylic alcohols

To a solution of benzylic alcohol (3.47 mmol, 1.00 eq.) in PhMe (3.5 mL) in a sealable microwave vial under N_2 was added *N,N*-diisopropylcarbamoyl chloride (681 mg, 4.16 mmol, 1.20 eq.) followed by Et_3N (0.63 mL, 4.51 mmol, 1.30 eq.). The vial was then sealed and heated under microwave irradiation at 150 °C for 2 h. The reaction mixture was then cooled to room temperature, filtered through a plug of SiO_2 with Et_2O as eluent, concentrated *in vacuo* and purified by bulb-to-bulb distillation under reduced pressure to give pure product.

Preparation of Secondary Benzoates



Synthesis of secondary alcohol from enantiopure epoxide

(*S*)-Propylene oxide (5.0 mL, 69 mmol, 1 eq., >99:1 *er*) was added dropwise to allyl magnesium bromide (1M in Et_2O , 103 mL, 103 mmol, 1.5 eq.) at -78 °C and the resulting mixture was stirred for 2 h at this temperature. The reaction mixture was then warmed to room temperature and sat. aq. NH_4Cl (100 mL) was added and the organic phase separated. The aqueous layer was washed with Et_2O (2×100 mL), the organic phases were combined, dried (MgSO_4) and concentrated *in vacuo* to afford (*S*)-hex-5-en-ol, >99:1 *er*.

Benzoate formation through Mitsunobu reaction

To a stirred solution of PPh_3 (11 mmol), secondary alcohol (10 mmol) and 2,4,6-triisopropylbenzoic acid (11.5 mmol) in THF (15 mL) at 0 °C (ice bath), was added DIAD (11 mmol) dropwise over 10 min. After stirring for 4 h at 0 °C (ice bath), the volatiles were removed *in vacuo*. The residue was dissolved in pentane (15 mL) and stirred in for 5 min. The white suspension was filtered and the filter cake washed with pentane (100 mL). The solvent was removed *in vacuo* and the residue purified by flash column chromatography to give pure benzoate.

Notes:

- Sometimes separation of the excess PPh_3 from the product by chromatography can be challenging. You can distil your benzoate or take it up in THF and add 30% aq. H_2O_2 to form PPh_3O , which is easier to separate.