

NMR Tube Washers from New Era



- Inert polymer and glass components
- Non-breakable polymer solvent tube
- All components replaceable
- Tube washer fits standard vacuum flask
- To wash, simply add solvent to reservoir
- Will air dry once solvent is depleted

Tube washers

Part No.	Description
NE-231-RC3	3mm NMR tube washer
NE-231-RC5	5mm NMR tube washer
NE-240	Vacuum flask, 500ml

Standard components

Part No.	Description
NE-231-B	Base only, PTFE
NE-231-R	Reservoir only, PTFE
NE-321-G	Glass body only

5mm components

Part No.	Description
NE-231-C5	Cap gasket, 5mm, PE
NE-231-N5	Nut, PPS, 1/8"
NE-231-F5	Ferrule, ETFE, 1/8"
NE-231-S5	Solvent capillary, PEEK, 1/8"

3mm components

Part No.	Description
NE231-C3	Cap gasket, 3mm, PE
NE-231-N3	Nut, PPS, 1/16"
NE-231-F3	Ferrule, ETFE, 1/16"
NE321-S3	Solvent capillary, PEEK, 1/16"